

ZNF329 Antibody

Catalog No: #35186

Package Size: #35186-1 50ul #35186-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF329 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF329 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF329.
Target Name	ZNF329
Other Names	Zinc finger protein 329; ZNF329;
Accession No.	Swiss-Prot: Q86UD4NCBI Gene ID: 79673
SDS-PAGE MW	61kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HuvEc cells, using ZNF329 antibody #35186.

Background

May be involved in transcriptional regulation.

Sugano S., Nat. Genet. 36:40-45(2004).

Lucas S.M., Nature 428:529-535(2004).

Venter J.C., Submitted (JUL-2005).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.