

ZNF695 Antibody

Catalog No: #35191

Package Size: #35191-1 50ul #35191-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF695 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF695 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human ZNF695.
Target Name	ZNF695
Other Names	Zinc finger protein 695; Zinc finger protein SBZF3; ZNF695;
Accession No.	Swiss-Prot: Q8IW36NCBI Gene ID: 57116
SDS-PAGE MW	20kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat/293 cells, using ZNF695 antibody #35191.

Background

May be involved in transcriptional regulation.

Gregory S.G., Nature 441:315-321(2006).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.