

## IL-7Ra CD127 Polyclonal Antibody Cy7 Conjugated

Catalog No: #C03360Cy7

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | IL-7Ra CD127 Polyclonal Antibody Cy7 Conjugated   |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Isotype               | IgG   |
| Purification          | Purified by Protein A.  |
| Applications          | IF(IHC-P) FC  |
| Species Reactivity    | Hu Ms   |
| Immunogen Description | KLH conjugated synthetic peptide aa 250-300 459 derived from human IL-7Ra   |
| Conjugates            | Cy7   |
| Target Name           | IL-7Ra CD127  |
| Other Names           | ILRA; CD127; IL7RA; CDW127; IL-7R-alpha; Interleukin-7 receptor subunit alpha; IL-7 receptor subunit alpha; IL-7R subunit alpha; IL-7RA; IL7R |
| Accession No.         | Swiss-Prot#:P16871NCBI Gene ID:3575   |
| Cell Localization     | Cytoplasm   |
| Concentration         | 1mg/ml  |
| Formulation           | Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.   |
| Storage               | Store at 4C for 12 months.  |

## Application Details

IF:1:50-200

Flow-Cyt 1 &amp;mu;g /test

## Background

Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.