

CD171 NCAM-L1 Polyclonal Antibody Cy3 Conjugated

Catalog No: #C03859Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD171 NCAM-L1 Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	FC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 1144-1200 1257 derived from human CD171
Conjugates	Cy3
Target Name	CD171 NCAM-L1
Other Names	S10; HSAS; MASA; MIC5; SPG1; CAML1; CD171; HSAS1; N-CAML1; NCAM-L1; N-CAM-L1; Neural cell adhesion molecule L1; L1CAM
Accession No.	Swiss-Prot#:P32004NCBI Gene ID:3897
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Flow-Cyt: 1µg /test

Background

Cell adhesion molecule with an important role in the development of the nervous system. Involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc. Binds to axonin on neurons.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.