

# Jak1 (Phospho-Tyr1034 + Tyr1035) Antibody Cy5 Conjugated

Catalog No: #C04568Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Jak1 (Phospho-Tyr1034 + Tyr1035) Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF(IHC-P)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic phosphopeptide derived from human Jak1 around the phosphorylation site of Tyr1034 + Tyr1035
Conjugates	Cy5
Target Name	Jak1 Tyr1034 + Tyr1035
Other Names	Tyrosine-protein kinase JAK1; Janus kinase 1; JAK-1; JAK1; JAK1A; JAK1B
Accession No.	NCBI Gene ID:3716
Concentration	1mg/ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

IF: 1:50-200

## Background

Tyrosine kinase of the non-receptor type, involved in the IFN-alpha beta gamma signal pathway. Kinase partner for the interleukin (IL)-2 receptor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.