

AMPK α 1/AMPK α 2(Phospho-Thr183/Thr172) Antibody

Catalog No: #11183



Package Size: #11183-1 50ul #11183-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	AMPK α 1/AMPK α 2(Phospho-Thr183/Thr172) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of AMPK α 1/AMPK α 2 only when phosphorylated at threonine 183 or 172.
Immunogen Type	Peptide-KLH
Immunogen Description	The antiserum was produced against synthesized peptide derived from human AMPK alpha around the phosphorylation site of Thr172.
Conjugates	Unconjugated
Target Name	AMPK α 1/AMPK α 2
Modification	Phospho
Other Names	AMPK; AMPK α 1; AMPK2
Accession No.	Swiss-Prot: Q13131/P54646; NCBI Gene ID: 5562/5563; NCBI mRNA: NM_006251.5/ NM_006252.3; NCBI Protein: NP_006242.5 /NP_006243.2
SDS-PAGE MW	63
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

WB 1:500-1:2000IHC 1:100-1:300IF 1:50-200;


Images



Western Blot analysis of various cells using Phospho-AMPK α 1/2 (T183/172) Polyclonal Antibody diluted at 1:500



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtained from antibody was pre-absorbed by immunogen peptide.



Immunofluorescence analysis of Hela cell. 1.AMPKα1/2 (phospho Thr183/172) Polyclonal Antibody(green) was diluted at 1:200(4° overnight). (red) was diluted at 1:200(4° overnight). 2. Goat Anti Rabbit Alexa Fluor 488 was diluted at 1:1000(room temperature, 50min). Goat Anti Mouse Alexa Fluor 594 was diluted at 1:1000(room temperature, 50min).

Background

Catalytic subunit of AMP-activated protein kinase (AMPK), an energy sensor protein kinase that plays a key role in regulating cellular energy metabolism. In response to reduction of intracellular ATP levels, AMPK activates energy-producing pathways and inhibits energy-consuming processes: inhibits protein, carbohydrate and lipid biosynthesis, as well as cell growth and proliferation. AMPK acts via direct phosphorylation of metabolic enzymes, and by longer-term effects via phosphorylation of transcription regulators. Also acts as a regulator of cellular polarity by remodeling the actin cytoskeleton; probably by indirectly activating myosin. Regulates lipid synthesis by phosphorylating and inactivating lipid metabolic enzymes such as ACACA, ACACB, GYS1, HMGCR and LIPE; regulates fatty acid and cholesterol synthesis by phosphorylating acetyl-CoA carboxylase (ACACA and ACACB) and hormone-sensitive lipase (LIPE) enzymes, respectively. Regulates insulin-signaling and glycolysis by phosphorylating IRS1, PFKFB2 and PFKFB3. AMPK stimulates glucose uptake in muscle by increasing the translocation of the glucose transporter SLC2A4/GLUT4 to the plasma membrane, possibly by mediating phosphorylation of TBC1D4/AS160. Regulates transcription and chromatin structure by phosphorylating transcription regulators involved in energy metabolism such as CRTC2/TORC2, FOXO3, histone H2B, HDAC5, MEF2C, MLXIPL/ChREBP, EP300, HNF4A, p53/TP53, SREBF1, SREBF2 and PPARGC1A. Acts as a key regulator of glucose homeostasis in liver by phosphorylating CRTC2/TORC2, leading to CRTC2/TORC2 sequestration in the cytoplasm. In response to stress, phosphorylates 'Ser-36' of histone H2B (H2BS36ph), leading to promote transcription. Acts as a key regulator of cell growth and proliferation by phosphorylating TSC2, RPTOR and ATG1: in response to nutrient limitation, negatively regulates the mTORC1 complex by phosphorylating RPTOR component of the mTORC1 complex and by phosphorylating and activating TSC2. In response to nutrient limitation, promotes autophagy by phosphorylating and activating ULK1. AMPK also acts as a regulator of circadian rhythm by mediating phosphorylation of CRY1, leading to destabilize it. May regulate the Wnt signaling pathway by phosphorylating CTNNB1, leading to stabilize it. Also has tau-protein kinase activity: in response to amyloid beta A4 protein (APP) exposure, activated by CAMKK2, leading to phosphorylation of MAPT/TAU; however the relevance of such data remains unclear in vivo. Also phosphorylates CFTR, EEF2K, KLC1, NOS3 and SLC12A1.

Published Papers

el at., Gemcitabine induces apoptosis and autophagy via the AMPK/mTOR signaling pathway in pancreatic cancer cells, In Biotechnol Appl Biochem. On 2018 Sep by Zhu J, Chen Y et al..PMID: 29575133, , (2018)

[PMID:29575133](#)

el at., Eupatilin ameliorates dextran sulphate sodium-induced colitis in mice partly through promoting AMPK activation. In Phytomedicine. On 2018 Jul 15 by Zhou K, Cheng R et al..PMID:30097122, , (2018)

[PMID:30097122](#)

el at., Concurrent Control of the Kaposi's Sarcoma-Associated Herpesvirus Life Cycle through Chromatin Modulation and Host Hedgehog Signaling: a New Prospect for the Therapeutic Potential of Lipoxin A4. In J Virol on 2020 Apr 16 by Asha K, Balfe N, et al.. PMID:32102879, , (2020)

[PMID:32102879](#)

el at., FSTL1 as a Potential Mediator of Exercise-Induced Cardioprotection in Post-Myocardial Infarction Rats. In Sci Rep on 2016 Aug 26 by Yue Xi, Da-Wei Gong et al.. PMID: 27561749, , (2016)

[PMID:27561749](#)

el at., Glutaredoxins concomitant with optimal ROS activate AMPK through S-glutathionylation to improve glucose metabolism in type 2 diabetes. In Free Radic Biol Med. On 2016 Dec by Dong K, Wu M et al.. PMID:27743883, , (2016)

[PMID:27743883](#)

el at., Neuronal Ca²⁺ sensor-1 contributes to stress tolerance in cardiomyocytes via activation of mitochondrial detoxification pathways. In J Mol Cell Cardiol on 2016 Oct by Tomoe Y Nakamura, Shu Nakao et al.. PMID:27555477, , (2016)

[PMID:27555477](#)

el at., MACC1 mediates acetylcholine-induced invasion and migration by human gastric cancer cells. In Oncotarget on 2016 Apr 5 by Ting Yang, Wanming He et al.. PMID: 26919111, , (2016)

[PMID:26919111](#)

el at., Hepatic AMPK activation in response to dynamic REDOX balance is a biomarker of exercise to improve blood glucose control. In Elife on 2022 Sep 26 by Meiling Wu, Anda Zhao, et al.. PMID:36155132, , (2022)

[PMID:36155132](#)

el at., Hepatic AMPK signaling dynamic activation in response to REDOX balance are sentinel biomarkers of exercise and antioxidant intervention to improve blood glucose control. In Elife on 2022 Sep 26 by Meiling Wu, Anda Zhao, et al.. PMID:36155132, , (2022)

[PMID:36155132](#)

el at., OR2H2 Activates CAMKK β AMPK β Autophagy Signaling Axis and Suppresses Senescence in VK2/E6E7 Cells. In Pharmaceuticals (Basel) On 2023 Sep by Ji Min Kim, Sina Dziobaka et al.. PMID:37765029, , (2023)

[PMID:37765029](#)

el at., The AMPK-dependent inhibition of autophagy plays a crucial role in protecting photoreceptor from photooxidative injury. In J Photochem Photobiol B On 2023 Aug by Yu-Lin Li, Tian-Zi Zhang et al.. PMID:37302163, , (2023)

[PMID:37302161](#)

el at., Docosahexaenoic Acid Coordinating with Sodium Selenite Promotes Paraptosis in Colorectal Cancer Cells by Disrupting the Redox Homeostasis and Activating the MAPK Pathway. In Nutrients on 2024 Jun 1 by Sheng Zhao, Yuzhou Meng, et al.. PMID:38892670, , (2024)

[PMID:38892670](#)

Anni Zheng;Anni Zheng;Anni Zheng;Sufei Song;Sufei Song;Sufei Song;Qiuling Xu;Qiuling Xu;Qiuling Xu;Tao Liu;Tao Liu;Tao Liu el at., Potential Mechanisms for Shenling Baizhu Powder in the Treatment of MAFLD Based on Network Pharmacology and Experimental Validation, , (2025)

[PMID:](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.