

Tau (Phospho-Ser202) Antibody

Catalog No: #11723

Package Size: #11723-1 50ul #11723-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Tau (Phospho-Ser202) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of Tau only when phosphorylated at serine 202.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around phosphorylation site of serine 202(P-G-S(p)-P-G) derived from Human Tau.
Target Name	Tau
Modification	Phospho
Other Names	MAPT; MTBT1; PHF-tau;
Accession No.	Swiss-Prot#: P10636; NCBI Gene#: 4137; NCBI Protein#: NP_058519.3.
SDS-PAGE MW	85kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500~1:1000

Images



Western blot analysis of extracts from JK cells (Lane 2), using Tau (Phospho-Ser202) Antibody #11723. The lane on the left is treated with antigen-specific peptide.

Background

Promotes microtubule assembly and stability, and might be involved in the establishment and maintenance of neuronal polarity. The C-terminus binds axonal microtubules while the N-terminus binds neural plasma membrane components, suggesting that tau functions as a linker protein between both. Axonal polarity is predetermined by tau localization (in the neuronal cell) in the domain of the cell body defined by the centrosome. The short isoforms allow plasticity of the cytoskeleton whereas the longer isoforms may preferentially play a role in its stabilization.

Goedert M., Proc. Natl. Acad. Sci. U.S.A. 85:4051-4055(1988).

Goedert M., EMBO J. 8:393-399(1989).

Lee G., Neuron 2:1615-1624(1989).

Published Papers

el at., Brain-derived neurotrophic factor protects against tau-related neurodegeneration of Alzheimer's disease.In Transl Psychiatry on 2016 Oct 4 by S-S Jiao , L-L Shen et al..PMID:27701410 , , (2016)

[PMID:27701410](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.