

# Gab 1 (Phospho-Tyr659) Polyclonal Antibody

Catalog No: #12359

Package Size: #12359-1 50ul #12359-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gab 1 (Phospho-Tyr659) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Phospho-Gab 1 (Y659) Polyclonal Antibody detects endogenous levels of Gab 1 protein only when phosphorylated at Y659.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Gab 1 around the phosphorylation site of Y659.
Conjugates	Unconjugated
Target Name	Gab 1
Modification	Phospho
Other Names	GAB1; GRB2-associated-binding protein 1; GRB2-associated binder 1; Growth factor receptor bound protein 2-associated protein 1
Accession No.	Swiss-Prot: Q13480NCBI Gene ID: 2549
Target Species	human
SDS-PAGE MW	76kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/5000

Not yet tested in other applications

## Images



Western Blot analysis of 3T3 cells using Phospho-Gab 1 (Y659) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.