

Androgen Receptor (Phospho-Tyr267) Antibody

Catalog No: #12557

Package Size: #12557-1 50ul #12557-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Androgen Receptor (Phospho-Tyr267) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Human
Specificity	Androgen Receptor (Phospho-Tyr267) Antibody detects endogenous levels of Androgen Receptor only when phosphorylated at Tyr267
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Androgen Receptor (Phospho-Tyr267)
Conjugates	Unconjugated
Target Name	Androgen Receptor
Modification	Phospho
Other Names	AR, Androgen receptor, Dihydrotestosterone receptor, DHTR, AIS, KD, NR3C4, HUMARA, SBMA, TFM, SMAX1, HYPSP1
Accession No.	Swiss-Prot#: P10275NCBI Gene ID: 367
Target Species	human
Calculated MW	110kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of 1. A549 cell and 2. A549 cell treated with 100ng/mL LPS for 30min, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.