

FoxM1 (Phospho-Thr600) Antibody

Catalog No: #13207

Package Size: #13207-1 50ul #13207-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FoxM1 (Phospho-Thr600) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho- and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	Phospho-FoxM1 (Phospho-Thr600) Antibody detects endogenous levels of FoxM1 only when phosphorylated at Thr600
Immunogen Description	A synthesized peptide derived from human FoxM1 around the phosphorylation site of Thr600
Conjugates	Unconjugated
Target Name	FoxM1
Other Names	FKHL16; Forkhead box M1; Forkhead box protein M1; forkhead like 16; Forkhead-related protein FKHL16; FOX M1; Foxm1; FOXM1_HUMAN; FOXM1B; Hepatocyte nuclear factor 3 forkhead homolog 11; HFH-11; HFH11; HNF-3/fork-head homolog 11; HNF3; INS1; M phase phosphoprotein 2; M-phase phosphoprotein 2; MPHOSPH2; MPM-2 reactive phosphoprotein 2; MPP2; PIG29; TGT3; Transcription factor Trident; Trident; WIN; Winged-helix factor from INS-1 cells; Winged-helix factor from INS1 cells;
Accession No.	Swiss-Prot#:Q08050NCBI Gene ID:2305
Calculated MW	83
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

WB: 1:500-1:2000

Published Papers

el at., The potent novel CDK4/6 inhibitor TQB3616 in hormone receptor positive breast cancer: preclinical characterization with in vitro and human tumor xenograft modelsInBreast Cancer?On2023 Dec 8byWenyu Hu?#?1,?Lei Wang et al..PMID:?38090281, , (2023)

[PMID:38090281](https://pubmed.ncbi.nlm.nih.gov/38090281/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.