

KIF1B (phospho-Ser1487) rabbit pAb

Catalog No: #13764

Package Size: #13764-1 50ul #13764-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KIF1B (phospho-Ser1487) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human KIF1B (phospho-Ser1487)
Immunogen Description	Synthesized phospho peptide around human KIF1B (Ser1487)
Other Names	Kinesin-like protein KIF1B (Klp)
Accession No.	Swiss Prot:O60333GeneID:23095
SDS-PAGE MW	200
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

This gene encodes a motor protein that transports mitochondria and synaptic vesicle precursors. Mutations in this gene cause Charcot-Marie-Tooth disease, type 2A1. [provided by RefSeq, Jul 2008].

Note: This product is for in vitro research use only and is not intended for use in humans or animals.