

YWHAZ Polyclonal Antibody

Catalog No: #21188



Package Size: #21188-1 50ul #21188-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	YWHAZ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB IHC
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of total 14-3-3 z protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.56~60 (R-S-S-W-R) derived from Human 14-3-3 zeta.
Conjugates	Unconjugated
Target Name	14-3-3z
Other Names	1433Z; 143Z; FAS; Factor activating exoenzyme S; KCIP-1
Accession No.	Swiss-Prot: P63104NCBI Protein: NP_001129171.1
Calculated MW	28kDa
SDS-PAGE MW	28kDa
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

WB 1:500-1:2000; IHC 1:100-1:300;

Images



Western blot analysis of extracts from HepG2 , HeLa and 293 cells using 14-3-3z(Ab-58) Antibody #21188.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using 14-3-3z(Ab-58) Antibody #21188(left) or the same antibody preincubated with blocking peptide(right).

Background

This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene.

Gu YM, et al. (2006) FEBS Lett ; 580(1): 305-310

Powell DW, et al. (2003) Mol Cell Biol; 23(15): 5376-5387

Mackintosh C. (2004) Biochem J; 381(Pt 2): 329-342.

Published Papers

Alba Ferre;Francois Chauvigne;Magdalena Gozdowska;Ewa Kulczykowska;Roderick Nigel Finn;Joan Cerda et al., Neurohypophysial and paracrine vasopressinergic signaling regulates aquaporin trafficking to hydrate marine teleost oocytes, , (2023)

[PMID:37635977](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.