

Tuftelin 1 antibody

Catalog No: #22036



Package Size: #22036 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tuftelin 1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 134 and 318 of Tuftelin 1
Target Name	Tuftelin 1
Accession No.	NCBI Gene ID: 7286NCBI mRNA#: BC002933NCBI Protein#: AAH02933
Concentration	0.8mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 42kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)
 A: JurKat
 10% SDS PAGE
 Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Colon ca, using Tuftelin 1 antibody at 1: 250 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.