

## IFIT3 antibody

Catalog No: #22040

Package Size: #22040 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	IFIT3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 230 and 470 of Human IFIT3
Target Name	IFIT3
Accession No.	NCBI Gene ID: 3437NCBI mRNA#: NM_001031683NCBI Protein#: NP_001026853
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


## Application Details

Predicted MW: 56kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

## Images




Sample (30 ug of whole cell lysate)

A: 293T

7.5% SDS PAGE

Primary antibody diluted at 1: 5000



Immunohistochemical analysis of paraffin-embedded Hela xenograft, using IFIT3 antibody at 1: 500 dilution.

## Published Papers

el at., NF- $\kappa$ B-Dependent IFIT3 Induction by HBx Promotes Hepatitis B Virus Replication. In *Front Microbiol* on 2019 Oct 11 by Xu F, Song H, et al..PMID:31681236, (2019)

[PMID:31681236](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.