

Interferon alpha-8 antibody

Catalog No: #22057

Package Size: #22057 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon alpha-8 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 174 of Human IFNA8
Target Name	Interferon alpha-8
Accession No.	NCBI Gene ID: 3445NCBI mRNA#: NM_002170NCBI Protein#: NP_002161
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 22kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Sample(30 ug whole cell lysate)
 A: Raji
 12% SDS PAGE
 Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded OVCA xenograft, using interferon alpha 8 antibody at 1: 500 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.