

Interferon alpha-6 antibody

Catalog No: #22059

Package Size: #22059 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Interferon alpha-6 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 125 and 189 of Interferon alpha-6
Target Name	Interferon alpha-6
Accession No.	NCBI Gene ID: 3443NCBI mRNA#: NM_021002NCBI Protein#: NP_066282
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 22kd

Western blotting: 1:100-1:1000

Images



Sample(30 ug of whole cell lysate)
A: H1299
12% SDS PAGE
Primary antibody diluted at 1: 200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.