

PHKA2 antibody

Catalog No: #22226

Package Size: #22226 100ul



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PHKA2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 1172 and 1235 of human PHKA2
Target Name	PHKA2
Other Names	PHK; PYK; XLG; PYKL; XLG2; GSD9A; MGC133071
Accession No.	NCBI Gene ID: 5256NCBI mRNA#: NM_000292NCBI Protein#: NP_000283
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 138kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)
A: Hep G2
5% SDS PAGE
Primary antibody diluted at 1: 500



Immunohistochemical analysis of paraffin-embedded Colon ca, using PHKA2 antibody at 1: 250 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.