

## SGTA antibody

Catalog No: #22393



Package Size: #22393 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	SGTA antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 6 and 313 of Human SGTA
Target Name	SGTA
Accession No.	NCBI Gene ID: 6449NCBI mRNA#: NM_003021NCBI Protein#: NP_003012
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 34kd

Western blotting: 1:500-1:3000

## Images



Sample (30 ug of whole cell lysate)

A: Hep G2

B: Molt-4

10% SDS PAGE

Primary antibody diluted at 1: 5000

## Background

This gene encodes a protein which is capable of interacting with the major nonstructural protein of parvovirus H-1 and 70-kDa heat shock cognate protein; however, its function is not known. Since this transcript is expressed ubiquitously in various tissues, this protein may serve a housekeeping function. [provided by RefSeq]

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.