

TCF19 antibody

Catalog No: #22706

Package Size: #22706 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	TCF19 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 296 and 359 of Human TCF19
Target Name	TCF19
Accession No.	NCBI Gene ID: 6941NCBI mRNA#: HSU25826NCBI Protein#: AAB60363
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 38kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)
A: 293T
10% SDS PAGE
Primary antibody diluted at 1: 1500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.