

TMEM59 antibody

Catalog No: #23096

Package Size: #23096 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	TMEM59 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 16 and 281 of Human TMEM59
Target Name	TMEM59
Accession No.	NCBI Gene ID: 9528NCBI mRNA#: NM_004872NCBI Protein#: NP_004863
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 36kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)

A: A431

12% SDS PAGE

Primary antibody diluted at 1: 500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.