

Pacific Blue anti-mouse CD45

Catalog No: #26846



Package Size: #26846-1 25T #26846-2 100T

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Pacific Blue anti-mouse CD45
Isotype	Rat IgG2b
Applications	FC
Species Reactivity	Mouse
Conjugates	Pacific Blue
Other Names	T200, Ly-5, LCA
Formulation	0.2%BSAε 0.1%Proclin 300η£·ι 'Eη • η'O ε ²ζI'AζA²oO PH7.2-7.4
Storage	4 ε'I'o'O ι• Ωε δΩ• ε-

Application Details

FC 5μl/ζ~'

Note: This product is for in vitro research use only and is not intended for use in humans or animals.