

APC anti-mouse CD206

Catalog No: #26850



Package Size: #26850-1 25T #26850-2 100T

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	APC anti-mouse CD206
Isotype	Rat IgG2a, k
Purification	FC
Species Reactivity	Mouse
Conjugates	APC
Formulation	0.2%BSAε 0.1%Proclin 300η£·ι 'Eη • η'O ε ²ζlAζA²οO PH7.2-7.4
Storage	4 εl'ο'O ι• Ωε δΩ• ε-

Application Details

FC 5μl/ζ~'

Note: This product is for in vitro research use only and is not intended for use in humans or animals.