

SMYD4 Monoclonal Antibody

Catalog No: #27215



Package Size: #27215-1 50ul #27215-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SMYD4 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7C8-B11-F10
Isotype	IgG1
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of SMYD4 and does not cross-react with related proteins.
Immunogen Type	Recombinant Protein
Immunogen Description	Purified recombinant human SMYD4 protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	SMYD4
Other Names	KIAA1936; SET and MYND domain containing 4; SET and MYND domain-containing protein 4; SMYD4; SMYD4_HUMAN; ZMYND21.
Accession No.	Uniprot: Q8IYR2 Gene ID: 114826
SDS-PAGE MW	89kd
Formulation	ascites
Storage	store at -20Λ C


Application Details

Western blotting: 1:200

Images



Western blot detection of SMYD4 antibody in 40ug HeLa,40ug MCF7 and 40ug MDA-MB-468 cell lysates using SMYD4 antibody (1:100 diluted).Predicted band size:89KDa.Observed band size:89KDa.



Western blot detection of SMYD4 antibody in HeLa, MCF7 and A431 cell lysates using SMYD4 antibody (1:200 diluted). Predicted band size: 89 kDa. Observed band size: 89 kDa.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.