

MESP1 Polyclonal Antibody

Catalog No: #27309



Package Size: #27309-1 50ul #27309-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MESP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human MESP1 (NP_061140.1).
Conjugates	Unconjugated
Other Names	MESP1; bHLHc5; mesoderm posterior protein 1
Accession No.	Swiss-Prot#:Q9BRJ9NCBI Gene ID:55897
Calculated MW	20kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:1000

Images



Western blot analysis of extracts of various cell lines, using MESP1 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.