

ITPRIP Polyclonal Antibody

Catalog No: #27566

Package Size: #27566-1 50ul #27566-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ITPRIP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse
Immunogen Description	Recombinant fusion protein of human ITPRIP (NP_203755.1).
Conjugates	Unconjugated
Other Names	ITPRIP; DANGER; KIAA1754; bA127L20; bA127L20.2; inositol 1,4,5-trisphosphate receptor-interacting protein
Accession No.	Swiss-Prot#:Q8IWB1NCBI Gene ID:85450
Calculated MW	62kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ITPRIP antibody.

Background

This gene encodes a membrane-associated protein that binds the inositol 1,4,5-trisphosphate receptor (ITPR). The encoded protein enhances the sensitivity of ITPR to intracellular calcium signaling. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.