

SPRR3 Polyclonal Antibody

Catalog No: #27599



Package Size: #27599-1 50ul #27599-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SPRR3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Mouse;Rat
Immunogen Description	Recombinant fusion protein of human SPRR3 (NP_005407.1).
Conjugates	Unconjugated
Other Names	SPRR3
Accession No.	Uniprot:Q9UBC9GeneID:6707
Calculated MW	26kDa
SDS-PAGE MW	26kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:1000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using SPRR3 antibody.



Immunofluorescence analysis of L929 cells using SPRR3 Rabbit pAb.

Background

Cross-linked envelope protein of keratinocytes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.