

VPS39 Polyclonal Antibody

Catalog No: #27889

Package Size: #27889-1 50ul #27889-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VPS39 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human VPS39 (NP_056104.2).
Conjugates	Unconjugated
Other Names	VPS39;TLP;VAM6;hVam6p
Accession No.	Uniprot:Q96JC1GeneID:23339
Calculated MW	102kDa
SDS-PAGE MW	102KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of Rat brain, using VPS39 antibody.

Background

This gene encodes a protein that may promote clustering and fusion of late endosomes and lysosomes. The protein may also act as an adaptor protein that modulates the transforming growth factor-beta response by coupling the transforming growth factor-beta receptor complex to the Smad pathway. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.