

MSL2 Polyclonal Antibody

Catalog No: #27900



Package Size: #27900-1 50ul #27900-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MSL2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human MSL2 (NP_060603.2).
Conjugates	Unconjugated
Other Names	MSL2; MSL-2; MSL2L1; RNF184; E3 ubiquitin-protein ligase MSL2
Accession No.	Swiss-Prot#:Q9HCI7NCBI Gene ID:55167
Calculated MW	62kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MSL2 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.