

Rat Anti-mouse F4/80/PE antibody

Catalog No: #28031



Package Size: #28031-1 50T #28031-2 100T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Anti-mouse F4/80/PE antibody
Host Species	Rat
Clonality	Monoclonal
Purification	Affinity purified by Protein G
Applications	FC
Species Reactivity	Mouse
Other Names	F4/80; Cell surface glycoprotein EMR1; Cell surface glycoprotein F4/80; DD7A5 7; Egf like module containing mucin like hormone receptor like 1; Egf like module containing mucin like hormone receptor like sequence 1; EGF like module receptor 1; EGF TM7; EGF-like module receptor 1; EGF-like module-containing mucin-like hormone receptor-like 1; EGFTM7; EMR 1; EMR-1; EMR1 hormone receptor; EMR1_HUMAN; Gpf480; Ly71; Lymphocyte antigen 71; TM7LN3.
Accession No.	Swiss-Prot#: Q61549 NCBI Gene ID:2015
Formulation	0.01M PBS, 0.5%BSA, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.

Application Details

Flow-Cyt=10ul/Test

Note: This product is for in vitro research use only and is not intended for use in humans or animals.