

CCDC70 Polyclonal Antibody

Catalog No: #28566

Package Size: #28566-1 50ul #28566-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CCDC70 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human CCDC70 (NP_112580.2).
Conjugates	Unconjugated
Other Names	CCDC70; coiled-coil domain containing 70
Accession No.	Swiss-Prot#:Q6NSX1NCBI Gene ID:83446
Calculated MW	29kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100

Images



Western blot analysis of extracts of rat liver, using CCDC70 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.