

ASZ1 Polyclonal Antibody

Catalog No: #28581



Package Size: #28581-1 50ul #28581-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ASZ1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human ASZ1 (NP_570124.1).
Conjugates	Unconjugated
Other Names	ASZ1; ALP1; ANKL1; C7orf7; CT1.19; GASZ; Orf3; ankyrin repeat, SAM and basic leucine zipper domain containing 1
Accession No.	Swiss-Prot#:Q8WWH4NCBI Gene ID:136991
Calculated MW	53kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:100


Images



Western blot analysis of extracts of various cell lines, using ASZ1 antibody.



Immunofluorescence analysis of rat testis cells using ASZ1 antibody.



Immunofluorescence analysis of mouse testis cells using ASZ1 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.