

NDNF Polyclonal Antibody

Catalog No: #28766



Package Size: #28766-1 50ul #28766-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NDNF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse;Rat
Immunogen Description	Recombinant fusion protein of human NDNF (NP_078850.3).
Conjugates	Unconjugated
Other Names	NDNF; C4orf31; NORD; protein NDNF
Accession No.	Swiss-Prot#:Q8TB73NCBI Gene ID:79625
Calculated MW	65kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NDNF antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.