

GIT1 Polyclonal Antibody

Catalog No: #28967



Package Size: #28967-1 50ul #28967-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GIT1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human GIT1 (NP_054749.2).
Conjugates	Unconjugated
Other Names	GIT1; GIT ArfGAP 1
Accession No.	Swiss-Prot#:Q9Y2X7NCBI Gene ID:28964
Calculated MW	100kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:200 - 1:2000IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using GIT1 antibody at 1:1000 dilution.




Immunofluorescence analysis of C6 cells using [KO] GIT1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of L929 cells using [KO] GIT1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U-2 OS cells using [KO] GIT1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts from normal (control) and GIT1 knockout (KO) HeLa cells, using GIT1 antibody at 1:3000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.