

ATOH8 Polyclonal Antibody

Catalog No: #29999



Package Size: #29999-1 50ul #29999-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ATOH8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse;Rat
Immunogen Description	A synthetic peptide of human ATOH8.
Conjugates	Unconjugated
Other Names	ATOH8; HATH6; bHLHa21; protein atonal homolog 8
Accession No.	Swiss-Prot#:Q96SQ7NCBI Gene ID:84913
Calculated MW	35kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ATOH8 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.