

# TRIP13 Polyclonal Antibody

Catalog No: #30024



Package Size: #30024-1 50ul #30024-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TRIP13 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human TRIP13 (NP_004228.1).
Conjugates	Unconjugated
Other Names	TRIP13; 16E1BP; pachytene checkpoint protein 2 homolog
Accession No.	Swiss-Prot#:Q15645NCBI Gene ID:9319
Calculated MW	48kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using TRIP13 antibody.

## Background

This gene encodes a protein that interacts with thyroid hormone receptors, also known as hormone-dependent transcription factors. The gene product interacts specifically with the ligand binding domain. This gene is one of several that may play a role in early-stage non-small cell lung cancer.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.