

MON1B Polyclonal Antibody

Catalog No: #30102



Package Size: #30102-1 50ul #30102-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MON1B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Rat
Immunogen Description	Recombinant fusion protein of human MON1B (NP_055755.1).
Conjugates	Unconjugated
Other Names	HSRG1, SAND2, SRG1
Accession No.	Swiss-Prot#:Q7L1V2NCBI Gene ID:22879
Calculated MW	59kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of U-87MG cells, using MON1B antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.