

MTCO1 Polyclonal Antibody

Catalog No: #30224



Package Size: #30224-1 50ul #30224-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MTCO1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant protein of human MTCO1.
Conjugates	Unconjugated
Other Names	COI; MTCO1; COX1
Accession No.	Swiss-Prot#:P00395NCBI Gene ID:4512
Calculated MW	40kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using MTCO1 antibody.



Immunofluorescence analysis of NIH-3T3 cells using MTCO1 Polyclonal Antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.