

pan Symmetric-DiMethyl-Arginine pAb

Catalog No: #30384



Package Size: #30384-1 50ul #30384-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	pan Symmetric-DiMethyl-Arginine pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	A synthetic peptide corresponding to a sequence containing symmetric dimethylated R.
Conjugates	Unconjugated
SDS-PAGE MW	22-70KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using pan Symmetric-DiMethyl-Arginine antibody.

Published Papers

Li Meng;Yucheng Jiang;Jiawen You;Yatao Chen;Shuyu Guo;Liming Chen;Junqing Ma et al., PRMT1-methylated MSX1 phase separates to control palate development., , (2025)
[PMID:39843447](https://pubmed.ncbi.nlm.nih.gov/39843447/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.