

DYNC1I1 Polyclonal Antibody

Catalog No: #30494



Package Size: #30494-1 50ul #30494-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DYNC1I1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human DYNC1I1 (NP_001129028.1).
Conjugates	Unconjugated
Other Names	DYNC1I1; DNCI1; DNCIC1; cytoplasmic dynein 1 intermediate chain 1
Accession No.	Swiss-Prot#:O14576NCBI Gene ID:1780
Calculated MW	73kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:200 - 1:2000IHC 1:20 - 1:200


Images



Western blot analysis of extracts of various cell lines, using DYNC1I1 antibody.



Immunohistochemistry of paraffin-embedded human stomach using DYNC111 antibody.



Immunohistochemistry of paraffin-embedded rat brain using DYNC111 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.