

NAA20 Polyclonal Antibody

Catalog No: #30824



Package Size: #30824-1 50ul #30824-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NAA20 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human NAA20 (NP_057184.1).
Conjugates	Unconjugated
Other Names	NAA20;NAT3;NAT3P;NAT5;NAT5P;dJ1002M8.1
Accession No.	Uniprot:P61599GenelD:51126
Calculated MW	20kDa
SDS-PAGE MW	20kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using NAA20 antibody.



Immunohistochemistry of paraffin-embedded rat kidney using NAA20 antibody.



Immunohistochemistry of paraffin-embedded human oophoroma using NAA20 antibody.

Background

NAT5 is a component of N-acetyltransferase complex B (NatB). Human NatB performs cotranslational N(alpha)-terminal acetylation of methionine residues when they are followed by asparagine (Starheim et al., 2008 [PubMed 18570629]).[supplied by OMIM, Apr 2009]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.