

ATG2B Polyclonal Antibody

Catalog No: #31572



Package Size: #31572-1 50ul #31572-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ATG2B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human ATG2B (NP_060506.5).
Conjugates	Unconjugated
Other Names	ATG2B; C14orf103; autophagy related 2B
Accession No.	Swiss-Prot#:Q96BY7NCBI Gene ID:55102
Calculated MW	270kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:1000

Images



Western blot analysis of extracts of various cell lines, using ATG2B antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.