

METTL14 Polyclonal Antibody

Catalog No: #31591



Package Size: #31591-1 50ul #31591-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	METTL14 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human METTL14 (NP_066012.1).
Conjugates	Unconjugated
Other Names	METTL14; hMETTL14; methyltransferase like 14
Accession No.	Swiss-Prot#:Q9HCE5NCBI Gene ID:57721
Calculated MW	65kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using METTL14 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.