

ZNF300 antibody

Catalog No: #31824

Package Size: #31824-1 50ul #31824-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF300 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Human
Immunogen Description	Fusion protein of human ZNF300
Conjugates	Unconjugated
Target Name	ZNF300
Accession No.	NCBI Protein#:BC117248
Concentration	0.4mg/ml
Formulation	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at -20°C/1 year


Application Details

IHC dilution:1:25-1:100

Images



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using (ZNF300 Antibody) at dilution 1/30.



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using (ZNF300 Antibody) at dilution 1/30.

Background

The protein encoded by this gene is a C2H2-type zinc finger DNA binding protein and likely transcriptional regulator. The function of this protein is not yet known. Three transcript variants encoding different isoforms have been found for this gene.

Published Papers

Zang Chencheng, Cai Min, Chen Qian, Yu Meng, Yan Jingjing, Guo Yilin, Wang Muqing, Wu Congyu, Gao Yuan, Qi Yun et al., Elephantopus scaber attenuates MC903-caused atopic dermatitis and decreases TNF- α /IFN- γ -induced chemokines by suppressing MKP-1-mediated AP-1 signaling β ¹, Journal of ethnopharmacology, (2025)

[PMID:40633861](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.