

RIMS3 antibody

Catalog No: #31983

Package Size: #31983-1 50ul #31983-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RIMS3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Fusion protein of human RIMS3
Conjugates	Unconjugated
Target Name	RIMS3
Other Names	NIM3; RIM3; RIM 3
Accession No.	NCBI Protein#:BC003103
Concentration	1mg/ml
Formulation	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at -20°C/1 year

Application Details

IHC dilution:1:50-1:200

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using (RIMS3 Antibody) at dilution 1/60.

Background

Regulates synaptic membrane exocytosis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.