

PTGES3 Antibody

Catalog No: #32773



Package Size: #32773-1 50ul #32773-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PTGES3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of total PTGES3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human PTGES3.
Conjugates	Unconjugated
Target Name	PTGES3
Other Names	P23; TEBP; cPGES;
Accession No.	Swiss-Prot:Q15185NCBI Gene ID:10728
SDS-PAGE MW	19KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PTGES3 antibody at 1:1000 dilution.

Background

Molecular chaperone that localizes to genomic response elements in a hormone-dependent manner and disrupts receptor-mediated transcriptional

activation, by promoting disassembly of transcriptional regulatory complexes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.