

JAK1 Antibody

Catalog No: #33721

Package Size: #33721-1 50ul #33721-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	JAK1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC IF
Species Reactivity	Human;Mouse
Specificity	The antibody detects endogenous levels of total JAK1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human JAK1.
Conjugates	Unconjugated
Target Name	JAK1
Other Names	EC 2.7.10.2; JAK-1; JAK1; JAK1A; Janus kinase 1
Accession No.	Swiss-Prot: P23458NCBI Gene ID: 3716
SDS-PAGE MW	130kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100


Immunofluorescence: 1:100~1:500

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using JAK1 antibody #33721.

Immunofluorescence analysis of HeLa cells, using JAK1 antibody #33721.



Background

Tyrosine kinase of the non-receptor type, involved in the IFN- α /beta/gamma signal pathway. Kinase partner for the interleukin (IL)-2 receptor.

Wilks A.F., Mol. Cell. Biol. 11:2057-2065(1991).

Lee S.-T., Oncogene 8:3403-3410(1993).

Lindholm C.K., Biochem. Biophys. Res. Commun. 296:929-936(2002).

Published Papers

Ighodaro Igbe;XiaoFei Shen;Wei Jiao;Zhe Qiang;Teng Deng;Sheng Li;WanLi Liu;HanWei Liu;GuoLin Zhang;Fei Wang et al., Dietary quercetin potentiates the antiproliferative effect of interferon- α in hepatocellular carcinoma cells through activation of JAK/STAT pathway signaling by inhibition of SHP2 phosphatase, , (2017)

[PMID:29371942](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.