

TR19L Antibody

Catalog No: #34020

Package Size: #34020-1 50ul #34020-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TR19L Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	The antibody detects endogenous levels of total TR19L protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human TR19L.
Conjugates	Unconjugated
Target Name	TR19L
Other Names	Tumor necrosis factor receptor superfamily member 19L [Precursor]; Receptor expressed in lymphoid tissues; RELT; TNFRSF19L;
Accession No.	Swiss-Prot: Q969Z4NCBI Gene ID: 84957
SDS-PAGE MW	46kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, COLO205 cells, 293 cells and HepG2 cells, using TR19L antibody #34020.

Background

Mediates activation of NF-kappa-B. May play a role in T-cell activation.

Sica G.L., Blood 97:2702-2707(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Cusick J.K., Biochem. Biophys. Res. Commun. 340:535-543(2006).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.