

# EBFP Rabbit Polyclonal Antibody

Catalog No: #38087

Package Size: #38087-1 50ul #38087-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                           |
|-----------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name          | EBFP Rabbit Polyclonal Antibody                                                                           |
| Host Species          | Rabbit                                                                                                    |
| Clonality             | Polyclonal                                                                                                |
| Purification          | Affinity purification using immunogen.                                                                    |
| Applications          | WB                                                                                                        |
| Specificity           | The EBFP Rabbit Polyclonal Antibody detects EBFP & EBFP Tag protein.                                      |
| Immunogen Description | Recombinant protein corresponding to EBFP-tag.                                                            |
| Conjugates            | Unconjugated                                                                                              |
| Target Name           | EBFP                                                                                                      |
| Other Names           | BFP; EBFP; Blue Fluorescent Protein; Enhanced Blue Fluorescent Protein                                    |
| Concentration         | 1.0mg/ml                                                                                                  |
| Formulation           | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage               | Store at -20°C                                                                                            |

## Application Details

WB 1:2000-1:5000

## Images



Western blot analysis of EBFP recombinant protein, using #38087 diluted at 1:5,000.

## Background

Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.