

AMBRA1 antibody

Catalog No: #38182



Package Size: #38182-1 50ul #38182-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	AMBRA1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous level of total AMBRA1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human Ambra1 (NP_060219.2).
Conjugates	Unconjugated
Target Name	AMBRA1
Other Names	AMBRA1;DCAF3;WDR94;FLJ20294
Accession No.	Uniprot:Q9C0C7GeneID:55626
SDS-PAGE MW	143kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using Ambra1 Antibody.



Immunohistochemistry of paraffin-embedded rat lung using AMBRA1 antibody.



Immunohistochemistry of paraffin-embedded mouse lung using AMBRA1 antibody.

Published Papers

Ye Liuqi;Lin Danlei;Zhang Wen;Chen Shiji;Zhen Yumiao;Akkouche Sara;Liang Xiaoxu;Chong Cheong-Meng;Zhong Hai-Jing; et al., AMBRA1 drives gastric cancer progression through regulation of tumor plasticity, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.