

IFIT2 Antibody

Catalog No: #39584



Package Size: #39584 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	IFIT2 Antibody
Purification	Affinity purified
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human IFIT2 protein 100-321aa (BC032839).
Target Name	IFIT2
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

WB: 1:200~1:1000
IHC: 1:100-1:500

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal liver showing cytoplasmic staining with anti-IFIT2 antibody at a dilution of 1/100.



Western Blot analysis of fetal liver and lung Lysate with anti-IFIT2 antibody

Published Papers

el at., Type III Interferon-Induced CBF β Inhibits HBV replication by hijacking HBx. In Cell Mol Immunol on 2019 Apr by Xu F, Song H, et al.. PMID: 29523836, , (2019)

[PMID:29523836](#)

el at., Prognostic Impact of Metabolism Reprogramming Markers Acetyl-CoA Synthetase 2 Phosphorylation and Ketohexokinase-A Expression In Non-Small-Cell Lung Carcinoma. In Front Oncol on 2019 Nov 5 by Yang X, Shao F, et al.. PMID:31750240, , (2019)

[PMID:31681236](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.